

3D VOLUMETRY & DIFFUSION TENSOR IMAGING POSTPROCESSING, APPLIED RESEARCH RESULTS

SZABÓ ADÁM GYÖRGY (1), GYEBNÁR GYULA(1), SIRÁLY ENIKŐ (2), FODOR ZSUZSANNA (2),
HIDÁSI ZOLTÁN (2), CSIBRI EVA (2), SALACZ PÁL (2), RUDÁS GÁBOR (1), CSUKLY GÁBOR (2)



(1) SEMMELWEIS UNIVERSITY MR RESEARCH CENTER
(2) SEMMELWEIS UNIVERSITY CLINIC OF PSYCHIATRY AND PSYCHOTHERAPY

MRKK DEMENTIA/MILD COGNITIVE IMPAIRMENT - POSSIBILITIES, OPTIONS

- Visual rating – sulci, ventricles, extracerebral liquorspaces – e.g. Scheltens-scale
- Planimetry – simple 2D measurements – the width of the Ill. ventricle or measuring the ventricle/brain ratio
- Volumetry – half automatized measurement to assess brain volumen data, cortical thickness, etc.
- Diffusion tensor imaging (DTI) – half automatized measurement to assess white matter data
- Resting state – half automatized measurement to assess functional connectivity impairment
- Automatized measurement to assess all above in the future...

MRKK DEMENTIA/MILD COGNITIVE IMPAIRMENT - POSSIBILITIES, OPTIONS

- Visual rating – sulci, ventricles, extracerebral liquorspaces – e.g. Scheltens-scale, Fazekas-scale
- Planimetry – simple 2D measurements – the width of the Ill. ventricle or measuring the ventricle/brain ratio
- Volumetry – half automatized measurement to assess brain volumen data, cortical thickness, etc.
- Diffusion tensor imaging (DTI) – half automatized measurement to assess white matter data
- Resting state – half automatized measurement to assess functional connectivity impairment
- Automatized measurement to assess all above in the future...

MRKK VOLUMETRY

Pros:

- It is (can be) detailed and accurate
- Good predictor in early stages as well


Cons:

- Slow (can be)
- Needs dedicated software and routine
- Database is needed and hard to interpret results on single subject


MRKK MOST OFTEN USED SOFTWARE


2 softwares are widely used in research

- SPM (voxel based morphometry toolbox)
 - Developed by UCL and Zürich universities
 - multimodal processing, analysis sw. (fMRI, PET, EEG)
 - Runs on windows, but the ui tedious and you better get familiar with Matlab
 - Faster than Freesurfer
- Freesurfer
 - Harvard developed
 - Needs Linux enviroment
 - You better get familiar with VMs if you don't dedicate a computer for it
 - Gives a wide type and amount of data
 - Slow (!)



VBM Analysis Pipeline Overview





VBM VS. FREESURFER


- VBM uses linear and non-linear deformation for intersubject analysis
- Freesurfer uses white matter geometry data (mesh) to inflate brain for intersubject analysis
- In essence VBM is Statistical Parametric Mapping of regional segmented tissue density or volume *
 - The exact interpretation of gray matter density or volume is complicated, and depends on the preprocessing steps used *
 - It is not interpretable as neuronal packing density or other cytoarchitectonic tissue properties *
 - The hope is that changes in these microscopic properties may lead to macro- or mesoscopic VBM-detectable differences

VBM VS. FREESURFER

- In Freesurfer one derives morphometric measures from geometric models of the cortical surface
- The yellow line is the surface boundary between cortical white matter and cortical gray matter known as the white surface; this represents the inner boundary of cortex.
- The red line is the boundary between the gray matter and dura and/or CSF; this is referred to as the pial surface.
- The cortex is modeled as a surface model which is a mesh of triangles as shown in the Figure. Each triangle is known as a face. The place where the corners of the triangles meet is called a vertex.

DIFFUSION WEIGHTED MRI

- Diffusion: Random (Brownian) motion of EC water molecules
- Hindering elements ↔ Tissue microstructure
 - CSF: \approx free, isotropic
 - Grey matter: hindered, isotropic
 - White matter: restricted, anisotropic
 - Early pathological changes on cellular level \leftrightarrow disruptions
- Dedicated MRI-sequence \rightarrow diffusion profile




DTI IN EARLY COGNITIVE DECLINE - GYULA GYEBNÁI

10

DTI – TENSOR MODEL

- Diffusion profile = ellipsoid (tensor)
 - 6 independent components
 - Gaussian process
- Tensor eigenvalues – shape and orientation
 - Color coding,
 - Tractography
- Scalar parameters:
 - Fractional anisotropy (**FA**: Shape – how stretched is the ellipsoid?)
 - Mean diffusivity (**MD**: Size – „volume“ of the ellipsoid)



DTI IN EARLY COGNITIVE DECLINE - GYULA GYEBNÁR

11

EXPLORE... – DTI

- Movement correction
- EPI correction (T1 registration)
- Susceptibility correction (induced distortion)
- Eddy current correction
- Tensor fitting
 - Outlier rejection



MILD COGNITIVE IMPAIRMENT=MCI

- Affects elderly most often
- Memory decline (noticed by patient or relatives)
- Mild decline in cognitive functions
- Self-care intact
- Dementia criteria are not met
- MCI most often precedes dementia
 - Conversion to dementia from MCI is about **10-15%**, compared to 1-4% in healthy elderly, which means they are **risk population** (Mitchell et al. Neurobiol Aging 2007)



MCI SUBTYPES

Amnesic MCI (aMCI): (memory impairment)		Non-amnesic MCI (naMCI): (executive, spatial-visual, language function impairment)	
Single-domain	Multi-domain	Single-domain	Multi-domain
Progression leads most often to Alzheimer type dementia		Progression leads most often to non Alzheimer type dementia (fronto-temporal dementia, Lewy-body dementia, Parkinson disease)	



AIMS

- To find a good clinical predictor
- To find potential central nervous structural differences that are helping clinicians to differentiate amnesic and non-amnesic subtype and to establish a proper diagnosis



METHODS

- Demographic and psychosocial status
- depression and anxiety tests (to exclude pseudodementia):
 - Geriatric Depression Scale (GDS)
 - Spielberger Trait- and Anxiety Scale (STAI)
- Cognitive function exploration:
 - Addenbrooke Cognitive Tests (AKV)
 - Mini Mental State Examination (MMSE) (exclude developed dementia)
 - Rey Auditory and Verbal Learning Test (RAVLT)
 - Trail Making Test A and B (TMT)
- MR scanning (3T Philips Achieva, 3D-T1, T2, Flair, DTI, T2-FFE, resting state)



SUBJECTS

62 subjects were involved

24 healthy control
Mean age: 65.5 years
(SD: 7.6)

18 naMCI
Mean age: 71 years
(SD: 7.3)

20 aMCI
Mean age: 71 years
(SD: 11.3)

- No significant differences in age, sex and education between the 3 groups ($p > 0.05$)
- **Excluding criteria:**
 - **dementia:** MMSE, age- and education specific cutoff values (min. limit: 21, max. limit: 27)
 - **depression:** GDS – no recent or actual depression episode (no sign. group diff. ($p > 0.1$))
 - trauma, epilepsy, mental retardation, psychotic episode or organic psychosyndrome




DIAGNOSIS

aMCI: Petersen criteria (Petersen et al. 2009)	naMCI
• Subjective memory decline	• Memory functions are intact
• Daily activity intact	
• Objective memory decline (Rey test I-5 or delayed remembrance, age specific cutoff values)	• Trail Making Test B or AKV sum (age, sex/education spec. cutoff values) • and VI/OM < 3.2 (AKV)
• Intact global mental functions (MMSE)	
• Dementia criteria not met	

MRK

MR EXAMINATION PARAMETERS

- 3D-T1 images were acquired on a 3T Philips MR
 - 3T Philips Achieva, sag. 3D-T1 (1 mm), axi. T2, FLAIR, DTI
- Processing: FreeSurfer
 - Motion correction
 - BET
 - Talairach transformation
 - Segmentation
 - Normalization
 - Topology correction
 - Cortical model fitting
- Group wise statistics: SAS



MRK

REGIONS OF INTEREST

- Hippocampus**
- Entorhinalis cortex**
- Gyrus parahippocampalis**
 - Gyrus fusiformis
 - Lobus parietalis superior
 - Lobus parietalis inferior
 - Gyrus cinguli posterior
- Isthmus cinguli**
- Precuneus**
 - Gyrus temporalis superior
 - Gyrus temporalis medius
 - Gyrus temporalis inferior
 - Gyrus supramarginalis
- Amygdala**

